

EXPERIMENT NO.: 5. A

DATE:

**AIM: TO FIND OUT BLOOD GROUP OF OWN BLOOD SAMPLE.**

**REQUIREMENT:** Blood group detection kit [antisera kit], glass slides, Permanent Marker Pen, 70% alcohol, Cotton swabs.



### PRINCIPLE

Compatibility between the blood groups of donor and recipient determines the success of a blood transfusion. The ABO and Rh blood groups are looked at while conducting the test. In a diagnostic lab, Monoclonal antibodies are available for A, B and Rh antigen. Monoclonal antibody against Antigen A (also called Anti-A), comes in a small bottles with droppers; the monoclonal suspension being **BLUE** in color. Anti-B comes in **YELLOW** colour. Anti-D (monoclonal antibody against Rh) is **COLORLESS**. All the colour codes are universal standards. When the monoclonal antibodies are added one by one to wells that contain the test sample (blood from patient), if the RBCs in that particular sample carry the corresponding Antigen, clumps can be observed in the corresponding wells. A drop of blood is left without adding any of the antibodies; it is used as a control in the experiment. The monoclonal antibody bottles should be stored in a refrigerator.

### THEORY


- The surface of the erythrocyte contains some glycoprotein and glycolipids that can act as antigen. These antigens are known as isoantigens or agglutinogens.
- Based on the presence or absence of various isoantigens blood is categorized in to different blood groups. The two most important ones are: ABO and the RhD antigen; they determine someone's blood type (A, B, AB and O, with +, – or Null denoting RhD status).

➤ **ABO blood groups:**

- The ABO blood group is based on two glycolipids isoantigens called A and B.

**Child blood type estimate table**

		Father's Blood Type			
		A	B	AB	O
Mother's Blood Type	<b>A</b>	A/O	A/B/AB/O	A/B/AB	A/O
	<b>B</b>	A/B/AB/O	B/O	A/B/AB	B/O
	<b>AB</b>	A/B/AB	A/B/AB	A/B/AB	A/B
	<b>O</b>	A/O	B/O	A/B	O
Blood Group		Antigens present		Antibodies present	
A		A Antigen		Anti-B	
B		B Antigen		Anti-A	
AB		A and B Antigens		No antibodies	
O		Neither Antigens		Anti-A and Anti-B	



➤ **Rh Blood groups:**

- Rh antigen first found in Rhesus monkey so it is known as Rh blood groups.
- People whose blood has Rh antigen is known as Rh positive (+).
- People whose blood does not have Rh antigen is known as Rh negative (-).
- According to Rh positive and Rh negative ABO blood group further divided into eight types:

1. A+ve blood group
2. B+ve blood group
3. AB+ve blood group
4. O+ve blood group
5. A-ve blood group
6. B-ve blood group
7. AB-ve blood group
8. O-ve blood group

- If Rh- person receives Rh+ blood, their immune system starts to make anti-Rh antibodies that remain in the blood and during the second transfusion the previously formed anti-Rh antibodies will cause hemolysis of donated blood and cause a severe reaction.

Example:


**BLOOD TYPE COMPATIBILITY:**

		Donor Blood Type							
		A+	A-	B+	B-	AB+	AB-	O+	O-
Recipient Blood Type	A+	√	√	X	X	X	X	√	√
	A-	X	√	X	X	X	X	X	√
	B+	X	X	√	√	X	X	√	√
	B-	X	X	X	√	X	X	X	√
	AB+	√	√	√	√	√	√	√	√
	AB-	X	√	X	√	X	√	X	√
	O+	X	X	X	X	X	X	√	√
	O-	X	X	X	X	X	X	X	√

**PROCEDURE**

1. Take neat and clean four glass slides, Mark A, B, and D on three slide top middle part respectively. Use fourth slide to mix the blood with antisera.
2. Sterilize the ring finger with spirit using cotton swab.
3. Prick the finger using lancet/pricking needle.
4. Discard the first drop of blood.
5. Put blood drop on Slide A, Slide B and Slide D at the center.
6. Take the Anti-A (blue) bottle, re-suspend the content and use the dropper to place a drop of the Anti-A in the Slide A. Place the bottle back in ice.
7. Take the Anti-B (yellow) bottle, re-suspend the content and use the dropper to place a drop of the Anti-B in the Slide B. Place the bottle back in ice.
8. Take the Anti-D (colorless) bottle, re-suspend the content and use the dropper to place a drop of the Anti-D in the Slide D. Place the bottle back in ice.
9. Take a fourth slide, mix the Slide A, Slide B and Slide D blood drops with Anti-A, Anti-B and Anti-D solution respectively using alternate corner of the fourth slide. So, it prevent the mixing of one slide solution with the other slide solution.
10. After mixing, wait for a while to observe clumps.




**DIAGRAM**



No clumping

Clumping

**OBSERVATION**

Blood drop + Anti A	Blood drop + Anti B	Blood drop + Anti RhD
		

Tick ( ✓ ) appropriate box.

Slide No. A:	Blood + Anti serum A	Clump [   ] / No Clump [   ]
Slide No. B:	Blood + Anti serum B	Clump [   ] / No Clump [   ]
Slide No. D:	Blood + Anti serum D (Rh Factor)	Clump [   ] / No Clump [   ]

**RESULT**

Blood group of own blood sample is \_\_\_\_\_

**SIGNATURE OF TEACHER**